



	Technical Data Sheet	
Product	Anti-Hu CD79a Purified	
Cat. Number/Size	11-869-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD79a	
Clone	ZL7.4	
Format	Purified	
Reactivity	Human	
Application	IHC-P, ICC, FC (QC tested), ELISA, IP	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody ZL7.4 interacts with extracellular domain of CD79a (Ig alpha), a 40-45 kDa subunit of B cell antigen-specific receptor (BCR) and its early developmental forms.	
Other names	BCR alpha, Ig-alpha, MB-1, IGA	
Workshop	HCDM listed	
Immunogen	IgM complex isolated from Daudi cells	
Entrez Gene ID	973	
Gene name	CD79A	
NCBI Full Gene Name	CD79a molecule	
UniProt ID	P11912	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.