



	Technical Data Sheet	
Product	Anti-Hu CD73 Purified	
Cat. Number/Size	11-675-C025	0.025 mg
	11-675-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD73	
Clone	AD2	
Format	Purified	
Reactivity	Human, Non-human primates	
Application	IHC-F, FC (QC tested), MC-CyTOF	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody AD2 recognizes CD73, a 70 kDa GPI-anchored 5'-nucleotidase expressed predominantly on the surface of T and B cell subsets, follicular dendritic cells and endothelial cells.	
Other names	ecto-5'-nucleotidase, NT5E, E5NT, CALJA, NTE	
Workshop	HLDA IV, HLDA V	
Immunogen	pre-B leukemia cells	
Entrez Gene ID	4907	
Gene name	NT5E	
NCBI Full Gene Name	5'-nucleotidase ecto	
UniProt ID	P21589	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.