



	Technical Data Sheet	
Product	Anti-Hu CX3CR1 Purified	
Cat. Number/Size	11-956-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CX3CR1	
Clone	2A9-1	
Format	Purified	
Reactivity	Human, Non-human primates	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml.	
Isotype	Rat IgG2b kappa	
Specificity	The rat monoclonal antibody 2A9-1 recognizes an extracellular epitope of CX3CR1, an important regulator of leukocyte trafficking and immune system balance.	
Other names	V28, CCRL1, GPR13, CMKDR1, GPRV28, CMKBRL1, fractalkine receptor	
Immunogen	human CX3CR1 transfectants	
Entrez Gene ID	1524	
Gene name	CX3CR1	
NCBI Full Gene Name	C-X3-C motif chemokine receptor 1	
UniProt ID	P49238	
Concentration	1 mg/ml	
Preparation	Purified by protein-G affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.