

Technical Data Sheet

Product	Anti-Ms CD5 PE	
Cat. Number/Size	1P-748-C025	0.025 mg
	1P-748-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD5
Clone	53-7.3
Format	PE
Reactivity	Mouse
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Excitation laser	blue (488 nm)
Isotype	Rat IgG2a kappa
Specificity	The rat monoclonal antibody 53-7.3 recognizes an extracellular epitope of CD5, a 67kDa single-chain transmembrane glycoprotein expressed on mature T lymphocytes, most of thymocytes and B-1 lymphocytes.
Other names	T1, LEU1
Immunogen	murine thymus or spleen cells
Entrez Gene ID	12507
Gene name	CD5
NCBI Full Gene Name	CD5 molecule
UniProt ID	P13379

Concentration	0.5 mg/ml
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.