

Technical Data Sheet

Product

Anti-BrdU Purified

Cat. Number/Size

11-682-C025

0.025 mg

11-682-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	BrdU
Clone	Bu20a
Format	Purified
Application	FC (QC tested), IHC(P), IHC(F)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG1
Specificity	The antibody Bu20a reacts specifically with BrdU incorporated into DNA during S-phase of a cell cycle. It is useful for detecting proliferating cells by flow cytometry or immunohistochemistry staining.
Other names	5-bromodeoxyuridine
Immunogen	Bromodeoxyuridine conjugated to BSA

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.