

Technical Data Sheet

Product

**Anti-Hu TCR alpha/beta APC**

Cat. Number/Size

**1A-607-T025**

**25 tests**

**1A-607-T100**

**100 tests**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	TCR alpha/beta
Clone	IP26
Format	APC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Excitation laser	red (633 nm)
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody IP26 recognizes a monomorphic extracellular determinant of TCR alpha/beta, the dominant subtype of T cell receptor expressed in human peripheral blood.
Other names	TCRA/B

---

Preparation	The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.

---

Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>
-----------------------	--

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.