

Technical Data Sheet

Product

Anti-Hu IgA Secretory Component Purified

Cat. Number/Size

11-112-C025 **0.025 mg**
11-112-C100 **0.1 mg**

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	IgA Secretory Component
Clone	SC-05
Format	Purified
Reactivity	Human
Application	IP, WB, IHC(P), IHC(F)
Application details	Immunohistochemistry (frozen sections): Examples of positive human tissues: secretory mucosa of gastrointestinal and respiratory tract, epithelia of salivary glands, endometrium, endocervix, kidney, prostate, mammary gland. Recommended dilution: 2-8 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody SC-05 reacts with 80 kDa human secretory component glycoprotein (both free and bound in secretory IgA). Recognized glycoprotein is specific membrane marker of glandular carcinomas.
Other names	Immunoglobulin A SC
Immunogen	Affinity-purified secretory component from human colostrum.

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.