



	Technical Data Sheet	
Product	<b>Anti-Hu CD127 Purified</b>	
Cat. Number/Size	<b>11-962-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	
Antigen	CD127	
Clone	rHIQB010	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The recombinant mouse monoclonal antibody rHIQB010 recognizes an extracellular epitope of human CD127, the alpha subunit of IL-7 receptor.	
Other names	IL7R, IMD104	
Entrez Gene ID	3575	
Gene name	IL7R	
NCBI Full Gene Name	interleukin 7 receptor	
UniProt ID	P16871	
Preparation	Purified by affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.