



	Technical Data Sheet	
Product	Anti-TPX2 Purified	
Cat. Number/Size	11-818-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	TPX2	
Clone	TPX2-01	
Format	Purified	
Reactivity	Rat, Mouse, Human	
Application	IP, WB, IHC-P, ICC, ELISA	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody TPX2-01 recognizes the epitope EPFVPKKEKKS (aa 636-646 in human) of TPX2, a microtubule-associated intracellular critical regulator of mitosis, which is overexpressed in many types of tumors and is a marker of worse prognosis in various cancers.	
Other names	DIL2, HCTP4, REPP86	
Immunogen	Human TPX2	
Entrez Gene ID	22974	
Gene name	TPX2	
NCBI Full Gene Name	TPX2 microtubule nucleation factor	
UniProt ID	Q9ULW0	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.