

Technical Data Sheet

Product	Anti-IFN gamma HRP	
Cat. Number/Size	1X-706-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	IFN gamma
Clone	4S.B3
Format	HRP
Reactivity	Non-human primates, Human
Application	ELISA (QC tested)
Application details	ELISA: This antibody is being used as detection antibody in combination with capture antibody NIB42.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 4S.B3 recognizes IFN-gamma, a 16-25 kDa cytokine produced by activated Th1 cells and NK cells. Binds both glycosylated and non-glycosylated protein.
Other names	Interferon gamma, IFN-gamma
Immunogen	Interferon gamma derived from human leukocytes
Entrez Gene ID	3458
Gene name	IFNG
NCBI Full Gene Name	interferon gamma
UniProt ID	P01579

Concentration	0.2 mg/ml
Preparation	Purified antibody is conjugated with activated horseradish peroxidase (HRP) of high specific enzyme activity under optimum conditions and unconjugated antibody and free HRP are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 0.01% thimerosal
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.