

Technical Data Sheet

Product	<b>Anti-Hu CD200 Purified</b>
Cat. Number/Size	<b>11-746-C100</b> <b>0.1 mg</b>
	<b>For Research Use Only.</b>
	<b>Not for use in diagnostic or therapeutic procedures.</b>

---

Antigen	CD200
Clone	OX-104
Format	Purified
Reactivity	Human
Application	FC (QC tested), IHC(P), IHC(F)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody OX-104 recognizes an extracellular epitope of CD200, a type-1 glycoprotein of the immunoglobulin superfamily, which is expressed in neurons, B and T cell subsets, keratinocytes, follicular dendritic cells, and ovarian cells.
Other names	MRC, MOX, OX-2
Workshop	HLDA VII: WS Code 70655
Immunogen	recombinant human CD200
Entrez Gene ID	4345
Gene name	CD200
NCBI Full Gene Name	CD200 molecule
UniProt ID	P41217

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.