



	Technical Data Sheet	
Product	Anti-Hu CD3 Biotin	
Cat. Number/Size	1B-514-C025	0.025 mg
	1B-514-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD3	
Clone	UCHT1	
Format	Biotin	
Reactivity	Non-human primates, Human	
Application	FC (QC tested), IP, IHC-F	
Application details	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 1-12 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The antibody UCHT1 recognizes an extracellular epitope on CD3 antigen of the TCR/CD3 complex on mature human T cells. The UCHT1 antibody reacts with the epsilon chain of the CD3 complex.	
Other names	CD3E, T3E, TCRE	
Workshop	HLDA I: WS Code T 3; HLDA III: WS Code T 126; HLDA III: WS Code T 471; HLDA VI: WS Code T 6T-CD3.1	
Immunogen	human thymocytes followed by Sezary T cells	
Entrez Gene ID	916	
Gene name	CD3E	
NCBI Full Gene Name	CD3 antigen, epsilon polypeptide	
UniProt ID	P07766	
Concentration	1 mg/ml	
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.