



Technical Data Sheet

Product	Anti-HLA-DR/DP Purified	
Cat. Number/Size	11-239-C025	0.025 mg
	11-239-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		
Antigen	HLA-DR/DP	
Clone	HL-38	
Format	Purified	
Reactivity	Pig, Human	
Application	IP, WB, IHC-P, FC (QC tested)	
Application details	Western blotting: Recommended dilution: 1 µg/ml Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG2a	
Specificity	The antibody HL-38 recognizes an extracellular common epitope on beta-chain of human HLA-DR and HLA-DP. DR and DP are the isotypes of human MHC Class II molecules expressed on antigen-presenting cells (APC; dendritic cells, B lymphocytes, monocytes, macrophages).	
Immunogen	Raji Burkitt's lymphoma cell line	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.