

Technical Data Sheet

Product

Anti-Hu CD62L Purified

Cat. Number/Size

11-449-C025

0.025 mg

11-449-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD62L
Clone	LT-TD180
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody LT-TD180 reacts with an extracellular epitope of CD62L (L-selectin), a 74-95 kDa single chain type I glycoprotein expressed on most peripheral blood B lymphocytes, T lymphocytes, monocytes and granulocytes; it is also present on a subset of NK cells and certain hematopoietic malignant cells.
Other names	L-selectin, LSEL, SELL, LAM-1, LECAM1, LEU8, LNHR, LYAM1, PLNHR
Immunogen	Peripheral blood leukocytes
Entrez Gene ID	6402
Gene name	SELL
NCBI Full Gene Name	selectin L
UniProt ID	P14151

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.