

Technical Data Sheet

Product	Anti-Hu CD229 PE	
Cat. Number/Size	1P-618-T025	25 tests
	1P-618-T100	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD229
Clone	HLY9.25
Format	PE
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody HLY9.25 (also known as HLY9.1.25) recognizes an extracellular epitope of CD229 / Ly9, a 100-120 kDa cell surface glycoprotein expressed on T and B cells.
Other names	SLAMF3, LY9
Workshop	HCDM listed
Immunogen	CD299-transfected 300.19 pre-B cell line
Entrez Gene ID	4063
Gene name	LY9
NCBI Full Gene Name	lymphocyte antigen 9
UniProt ID	Q9HBG7

Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.