

Technical Data Sheet

Product	Anti-Hu CD158z APC	
Cat. Number/Size	1A-945-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD158z
Clone	CH21
Format	APC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 3-12 µg/ml
Excitation laser	red (633 nm)
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody CH21 recognizes an extracellular epitope of human CD158z (KIR3DL3), a transmembrane glycoprotein of killer cell inhibitory receptor family.
Other names	KIR3DL3
Immunogen	human CD158z transfectants
Entrez Gene ID	115653
Gene name	KIR3DL3
NCBI Full Gene Name	killer cell immunoglobulin like receptor, three Ig
UniProt ID	Q8N743

Concentration	0.1 mg/ml
Preparation	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.