



	Technical Data Sheet	
Product	Anti-Ms CD45 PerCP	
Cat. Number/Size	PC-553-C025	0.025 mg
	PC-553-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD45	
Clone	EM-05	
Format	PerCP	
Reactivity	Mouse	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 2-4 µg/ml; positive control: murine peripheral blood leukocytes.	
Excitation laser	blue (488 nm)	
Isotype	Rat IgG	
Specificity	The rat monoclonal antibody EM-05 reacts with an extracellular epitope of murine CD45 antigen (Leukocyte Common Antigen), a single chain type I transmembrane protein expressed at high level on cells of hematopoietic origin, except erythrocytes and platelets.	
Other names	LCA, T200, LY5, GP180, TPC	
Immunogen	Murine peripheral blood leukocytes	
Entrez Gene ID	19264	
Gene name	Ptprc	
NCBI Full Gene Name	protein tyrosine phosphatase, receptor type, C	
UniProt ID	Q8C6Q7	
Concentration	0.5 mg/ml	
Preparation	Purified antibody is conjugated with activated Peridinin-Chlorophyll Protein (PerCP) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.