

Technical Data Sheet

Product	Anti-Rt CD4 Purified	
Cat. Number/Size	11-679-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD4
Clone	OX-35
Format	Purified
Reactivity	Rat
Application	FC (QC tested), IHC, ICC
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2a kappa
Specificity	The mouse monoclonal antibody OX-35 reacts with an extracellular epitope of rat CD4 transmembrane glycoprotein (55 kDa).
Other names	T4/Leu-3, L3T4
Immunogen	MLR generated rat Th cells
Entrez Gene ID	24932
Gene name	CD4
NCBI Full Gene Name	CD4 molecule
UniProt ID	P05540

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.