



	Technical Data Sheet	
Product	Anti-Hu CD370 Purified	
Cat. Number/Size	11-787-C025	0.025 mg
	11-787-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD370	
Clone	8F9	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG2a	
Specificity	The mouse monoclonal antibody 8F9 recognizes an extracellular epitope of CD370 / CLEC9A (DNGR1), a type II transmembrane protein functioning as an endocytic receptor on BDCA31+ dendritic cells and on a subset of CD14+ CD16- monocytes.	
Other names	CLEC9A, DNGR1, DNGR-1, UNQ9341	
Workshop	HLDA X	
Immunogen	RBL-2H3 cells expressing human CLEC9A fused to an HA epitope	
Entrez Gene ID	283420	
Gene name	CLEC9A	
NCBI Full Gene Name	C-type lectin domain containing 9A	
UniProt ID	Q6UXN8	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.