



	Technical Data Sheet	
Product	Anti-Hu CD245 Purified	
Cat. Number/Size	11-941-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD245	
Clone	DY12	
Format	Purified	
Reactivity	Human	
Application	WB, FC (QC tested), IHC, IP	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody DY12 recognizes an extracellular epitope of CD245, a transmembrane protein expressed mainly by peripheral blood lymphocytes, NK cells, monocytes, granulocytes and platelets.	
Other names	p220/240, DY12, DY35	
Workshop	HLDA VI, HLDA VII	
Immunogen	human NK-type YT2C2 cells	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.