



	Technical Data Sheet	
Product	Anti-LIME1 Purified	
Cat. Number/Size	11-410-C025	0.025 mg
	11-410-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	LIME1	
Clone	PAb (410)	
Format	Purified	
Reactivity	Human	
Negative species	Mouse	
Application	WB, IHC-P	
Isotype	Rabbit polyclonal	
Specificity	The polyclonal antibody reacts with the cytoplasmic domain of LIME, a 30 kDa Lck-interacting transmembrane adaptor expressed by T cells.	
Other names	LITA1	
Immunogen	Bacterially expressed intracellular fragment corresponding to aa 141-295 of human LIME.	
Entrez Gene ID	54923	
Gene name	LIME1	
NCBI Full Gene Name	Lck interacting transmembrane adaptor 1	
UniProt ID	Q9H400	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.