

Technical Data Sheet

Product	Anti-DAXX Purified	
Cat. Number/Size	11-351-C025	0.025 mg
	11-351-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	DAXX
Clone	DAXX-01
Format	Purified
Reactivity	Human
Application	WB, ICC
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody DAXX-01 recognizes Daxx, a cytoplasmic death domain containing protein mainly expressed in fetal and adult human and mouse tissue.
Other names	DAP6, EAP1, BING2
Immunogen	Recombinant C-terminal part (aa 558-740) of human Daxx.
Entrez Gene ID	1616
Gene name	DAXX
NCBI Full Gene Name	death domain associated protein
UniProt ID	Q9UER7

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.