

Technical Data Sheet

Product	Anti-Ms CD16/CD32 Purified	
Cat. Number/Size	11-141-C025	0.025 mg
	11-141-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD16/CD32
Clone	93
Format	Purified
Reactivity	Mouse
Application	FC (QC tested), IP
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Rat IgG2a lambda
Specificity	The rat monoclonal antibody 93 recognizes a common extracellular epitope of murine CD16 (Fcγ ₃ R1) and CD32 (Fcγ ₂ R1), the low affinity receptors for IgG.
Other names	CD16, CD32, Fcγ ₃ R1, Fcγ ₂ R1, FCGR3, FCGR2
Immunogen	Murine pre-B cells

Concentration	1 mg/ml
Preparation	Purified by protein-G affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.