



Technical Data Sheet

Product	<b>Anti-Phosphotyrosine Purified</b>	
Cat. Number/Size	<b>11-263-C025</b>	<b>0.025 mg</b>
	<b>11-263-C100</b>	<b>0.1 mg</b>
<b>For Research Use Only.</b>		
<b>Not for use in diagnostic or therapeutic procedures.</b>		
Antigen	Phosphotyrosine	
Clone	P-Tyr-01	
Format	Purified	
Reactivity	Species independent	
Application	WB, ICC, FC (QC tested)	
Application details	Western blotting: Recommended dilution: 1-2 µg/ml. Flow cytometry: Recommended dilution: 2-5 µg/ml. Intracellular staining.	
Isotype	Mouse IgG1	
Specificity	The antibody P-Tyr-01 detects tyrosine phosphorylation in activated cells.	
Other names	PY	
Immunogen	Phosphotyrosine conjugated to bovine serum albumin.	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.