

Technical Data Sheet

Product	Anti-HLA Class I Purified	
Cat. Number/Size	11-243-C025	0.025 mg
	11-243-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	HLA Class I
Clone	MEM-81
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Isotype	Mouse IgG1
Specificity	The antibody MEM-81 reacts with an extracellular epitope on human classical MHC Class I molecules in native cell-surface forms. MHC Class I molecules (MHC Class Ia) are expressed on the surface of human nucleated cell types.
Immunogen	Leukocytes of a patient suffering with LGL type leukemia.

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.