

Technical Data Sheet

Product	<b>Anti-Ms CD45R Purified</b>	
Cat. Number/Size	<b>11-715-C025</b>	<b>0.025 mg</b>
	<b>11-715-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD45R
Clone	RA3-6B2
Format	Purified
Reactivity	Cat, Mouse, Human
Application	FC (QC tested), IP, IHC-P, IHC-F, ICC
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Rat IgG2a kappa
Specificity	The rat monoclonal antibody RA3-6B2 recognizes an extracellular epitope on CD45R, which is expressed at all developmental stages of B cells, including activated B cells, but also on subsets of NK and T cells. T cells detected by this antibody are supposed to be in activated state.
Other names	PTPRCR, T200R
Immunogen	Abelson murine leukemia virus-induced pre-B tumor cells

---

Concentration	1 mg/ml
Preparation	Purified by protein-G affinity chromatography followed by size-exclusion chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.