

Technical Data Sheet

Product	Anti-HLA-G FITC	
Cat. Number/Size	1F-438-C025	0.025 mg
	1F-438-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	HLA-G
Clone	01G
Format	FITC
Reactivity	Human
Negative species	Mouse, Rat
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 5 µg/ml.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The antibody 01G recognizes membrane-bound form of HLA-G (full-length HLA-G1), but not soluble forms. HLA-G belongs to the MHC Class I molecules (MHC Class Ib; nonclassical) and it is expressed on the surface of trophoblast cells.
Immunogen	HLA-B27 transgenic mice were immunized with H-2 identical murine cells transfected with and expressing genes encoding HLA-G and human beta2-microglobulin.
Entrez Gene ID	3135
Gene name	HLA-G
NCBI Full Gene Name	major histocompatibility complex, class I, G
UniProt ID	P17693

Concentration	1 mg/ml
Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.