

Technical Data Sheet

Product	<b>Anti-TNF alpha Purified</b>	
Cat. Number/Size	<b>11-688-C025</b>	<b>0.025 mg</b>
	<b>11-688-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	TNF alpha
Clone	MAb11
Format	Purified
Reactivity	Pig, Human, Non-human primates
Application	FC (QC tested), IHC(F), ICC, ELISA
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody MAb11 recognizes human 17-26 kDa cytokine TNF-alpha (tumor necrosis factor alpha).
Other names	TNFA, DIF, TNLG1F
Immunogen	Recombinant human TNF-alpha
Entrez Gene ID	7124
Gene name	TNF
NCBI Full Gene Name	tumor necrosis factor
UniProt ID	P01375

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.