



|                       |                                                                                                                                                                                                                                                                          |                  |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                       | Technical Data Sheet                                                                                                                                                                                                                                                     |                  |
| Product               | <b>Anti-Hu CD64 APC-Cy™7</b>                                                                                                                                                                                                                                             |                  |
| Cat. Number/Size      | <b>T4-644-T025</b>                                                                                                                                                                                                                                                       | <b>25 tests</b>  |
|                       | <b>T4-644-T100</b>                                                                                                                                                                                                                                                       | <b>100 tests</b> |
|                       | <b>For Research Use Only.</b>                                                                                                                                                                                                                                            |                  |
|                       | <b>Not for use in diagnostic or therapeutic procedures.</b>                                                                                                                                                                                                              |                  |
| Antigen               | CD64                                                                                                                                                                                                                                                                     |                  |
| Clone                 | 10.1                                                                                                                                                                                                                                                                     |                  |
| Format                | APC-Cy™7                                                                                                                                                                                                                                                                 |                  |
| Reactivity            | Human, Non-human primates                                                                                                                                                                                                                                                |                  |
| Application           | FC (QC tested)                                                                                                                                                                                                                                                           |                  |
| Application details   | Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.                                               |                  |
| Excitation laser      | red (633 nm)                                                                                                                                                                                                                                                             |                  |
| Isotype               | Mouse IgG1 kappa                                                                                                                                                                                                                                                         |                  |
| Specificity           | The mouse monoclonal antibody 10.1 recognizes an extracellular epitope on CD64/FcγRI, a 72 kDa single chain type I glycoprotein, that is expressed on monocytes/macrophages, dendritic cells, and activated granulocytes. The epitope is sensitive to formalin fixation. |                  |
| Other names           | FcRI, IGFR1, FcγR1A                                                                                                                                                                                                                                                      |                  |
| Workshop              | HLDA III: WS Code M-250                                                                                                                                                                                                                                                  |                  |
| Immunogen             | Rheumatoid synovial fluid cells and fibronectin purified human monocytes                                                                                                                                                                                                 |                  |
| Entrez Gene ID        | 2209                                                                                                                                                                                                                                                                     |                  |
| Gene name             | FCGR1A                                                                                                                                                                                                                                                                   |                  |
| NCBI Full Gene Name   | Fc fragment of IgG receptor Ia                                                                                                                                                                                                                                           |                  |
| UniProt ID            | P12314                                                                                                                                                                                                                                                                   |                  |
| Preparation           | Purified antibody is conjugated with activated tandem dye of allophycocyanin-cyanine 7 (APC-Cy™7) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.                                                 |                  |
| Formulation           | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                                                                                                                                                                                  |                  |
| Storage and handling  | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.                                                                                                                                                                                                 |                  |
| Images and References | www.exbio.cz                                                                                                                                                                                                                                                             |                  |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.

Cy™ and CyDye™ are registered trademarks of Cytiva.