

Technical Data Sheet

Product	<b>Anti-GABA(B) R2 Purified</b>	
Cat. Number/Size	<b>11-551-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	GABA(B) R2
Clone	PAb (551)
Format	Purified
Reactivity	Rat, Mouse
Application	WB, ICC
Isotype	Rabbit polyclonal
Specificity	The polyclonal antibody recognizes C-terminus of gamma-aminobutyric acid (GABA) B receptor 2 (recognized epitope: the last 23 aa). GB2 apparent MW ~120 kDa.
Other names	GABABR2, GABBR2
Immunogen	Synthetic peptide (coupled with THG) derived from the last 23 aa of mouse GABA B receptor 2. 100% homology with human GB2.

---

Concentration	1 mg/ml
Preparation	Purified from rabbit serum by affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.