



Technical Data Sheet

Product	Anti-FTT0571 (F. tularensis ssp. tularensis) Purified	
Cat. Number/Size	11-181-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	FTT0571 (F. tularensis ssp. tularensis)
Clone	PAb (181)
Format	Purified
Reactivity	Francisella tularensis ssp. tularensis
Application	WB, ELISA
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody to FTT0571 reacts with Francisella tularensis subsp. tularensis antigen FTT0571 (FTT_0571). The antibody does not react with Francisella tularensis subsp. holarctica homolog.
Other names	FTT_0571
Immunogen	Recombinant FTT0571 protein produced in E. coli
Entrez Gene ID	3191185
Gene name	FTT_0571
NCBI Full Gene Name	FTT0571 protein
UniProt ID	Q5NH96
Concentration	1 mg/ml
Preparation	Purified by ligand affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.
Images and References	www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.