

Technical Data Sheet

Product	Anti-Hu CD11b Biotin	
Cat. Number/Size	1B-681-C025	0.025 mg
	1B-681-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD11b
Clone	ICRF44
Format	Biotin
Reactivity	Human, Non-human primates
Application	FC (QC tested), IP, IHC(F), ICC
Application details	Flow cytometry: Recommended dilution: 1-12 µg/ml.
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody ICRF44 recognizes an extracellular epitope of CD11b (Mac-1alpha), a 165-170 kDa type 1 transmembrane protein mainly expressed on monocytes, granulocytes and NK-cells.
Other names	Mac-1, Integrin alpha M, ITGAM, CR3A, MO1A, MAC1A
Workshop	HLDA IV: WS Code M047
Immunogen	Rheumatoid synovial cells and human monocytes.
Entrez Gene ID	3684
Gene name	ITGAM
NCBI Full Gene Name	integrin subunit alpha M
UniProt ID	P11215

Concentration	1 mg/ml
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.