



	Technical Data Sheet	
Product	Anti-Hu CD229 Purified	
Cat. Number/Size	11-618-C025	0.025 mg
	11-618-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD229	
Clone	HLy9.25	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB, ICC	
Application details	Flow cytometry: Recommended dilution: 6 µg/ml; positive control: peripheral blood lymphocytes.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody HLy9.25 (also known as HLy9.1.25) recognizes an extracellular epitope of CD229 / Ly9, a 100-120 kDa cell surface glycoprotein expressed on T and B cells.	
Other names	SLAMF3, LY9	
Workshop	HCDM listed	
Immunogen	CD299-transfected 300.19 pre-B cell line	
Entrez Gene ID	4063	
Gene name	LY9	
NCBI Full Gene Name	lymphocyte antigen 9	
UniProt ID	Q9HBG7	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.