

Technical Data Sheet

Product	<b>Anti-Hu CD231 Purified</b>	
Cat. Number/Size	<b>11-200-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD231
Clone	B2D
Format	Purified
Reactivity	Human
Application	FC (QC tested), IP, WB, IHC
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody B2D recognizes CD231 (TALLA-1, tetraspanin 7), a transmembrane glycoprotein expressed in neuronal tissue and T-ALL.
Other names	TSPAN7, TALLA-1, MXS1, MRX58, TM4SF2
Immunogen	Human T-ALL cell line THP-6

---

Concentration	1 mg/ml
Preparation	Purified from cell culture supernatant by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.