

Technical Data Sheet

Product	Anti-mCherry purified
Cat. Number/Size	11-918-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	mCherry
Clone	PAb (918)
Format	Purified
Application	WB, ICC, FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml, extracellular staining or intracellular staining - depending on expression. Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody PAb (918) recognizes a red fluorescent protein tag mCherry.
Immunogen	mCherry protein from Anaplasma marginale
Gene name	mCherry
UniProt ID	X5DSL3

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.