

Technical Data Sheet

Product	Anti-MAP2 Purified	
Cat. Number/Size	11-473-C025	0.025 mg
	11-473-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	MAP2
Clone	MT-08
Format	Purified
Reactivity	Pig, Human, Mouse
Application	IP, WB, IHC(P), IHC(F), ICC, ELISA
Isotype	Mouse IgG1
Specificity	The antibody MT-08 recognizes an epitope (aa 1375-1395) located in central domain of molecule Microtubule Associated Protein 2ab (MAP2ab), an intracellular antigen.
Other names	Microtubule-associated protein 2, MAP-2
Immunogen	Microtubule protein (bovine brain) enriched for kinesin
Entrez Gene ID	4133
Gene name	MAP2
NCBI Full Gene Name	microtubule associated protein 2
UniProt ID	P11137

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.