

Technical Data Sheet

Product	Anti-Hu CD352 Purified	
Cat. Number/Size	11-865-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD352
Clone	hsF6.4.20
Format	Purified
Reactivity	Human
Application	ICC, IP, FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-12 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody hsF6.4.20 recognizes an extracellular epitope of CD352, a transmembrane glycoprotein expressed on NK, T and B cells.
Other names	SLAMF6, KALI, NTBA, Ly108, SF2000
Immunogen	CD352 transfectants
Entrez Gene ID	114836
Gene name	SLAMF6
NCBI Full Gene Name	SLAM family member 6
UniProt ID	Q96DU3

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.