



	Technical Data Sheet	
Product	Anti-Hu IL-17A Purified	
Cat. Number/Size	11-937-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	IL-17A	
Clone	9F9	
Format	Purified	
Reactivity	Human	
Application	IP, ELISA, FC (QC tested), ICC	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml. Intracellular staining.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 9F9 recognizes human interleukin 17A (IL-17A; secreted or intracellular).	
Other names	Interleukin 17	
Immunogen	mammalian-derived human IL-17-IgG fusion protein, boost with recombinant human IL-17A	
Entrez Gene ID	3605	
Gene name	IL17A	
NCBI Full Gene Name	interleukin 17A	
UniProt ID	Q16552	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.