

Technical Data Sheet

Product	Anti-Hu CD328 Purified	
Cat. Number/Size	11-799-C025	0.025 mg
	11-799-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD328
Clone	6-434
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody 6-434 recognizes an extracellular epitope of CD328 (Siglec-7), a 75 kDa transmembrane glycoprotein expressed mainly on NK cells, dendritic cells and monocytes.
Other names	SIGLEC7, p75/AIRM1, D-siglec
Workshop	HLDA VIII: WS Code ARN-80652
Immunogen	human dendritic cells
Entrez Gene ID	27036
Gene name	SIGLEC7
NCBI Full Gene Name	sialic acid binding Ig like lectin 7
UniProt ID	Q9Y286

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.