



|                       |                                                                                                                                                                                                                                                                      |                  |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                       | Technical Data Sheet                                                                                                                                                                                                                                                 |                  |
| Product               | <b>Anti-Hu CD1a PE</b>                                                                                                                                                                                                                                               |                  |
| Cat. Number/Size      | <b>1P-174-T025</b>                                                                                                                                                                                                                                                   | <b>25 tests</b>  |
|                       | <b>1P-174-T100</b>                                                                                                                                                                                                                                                   | <b>100 tests</b> |
|                       | <b>For Research Use Only.</b>                                                                                                                                                                                                                                        |                  |
|                       | <b>Not for use in diagnostic or therapeutic procedures.</b>                                                                                                                                                                                                          |                  |
| Antigen               | CD1a                                                                                                                                                                                                                                                                 |                  |
| Clone                 | SK9                                                                                                                                                                                                                                                                  |                  |
| Format                | PE                                                                                                                                                                                                                                                                   |                  |
| Reactivity            | Human                                                                                                                                                                                                                                                                |                  |
| Application           | FC (QC tested)                                                                                                                                                                                                                                                       |                  |
| Application details   | Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.                                            |                  |
| Excitation laser      | blue (488 nm)                                                                                                                                                                                                                                                        |                  |
| Isotype               | Mouse IgG2b kappa                                                                                                                                                                                                                                                    |                  |
| Specificity           | The mouse monoclonal antibody SK9 recognizes an extracellular epitope of CD1a (T6), a 49 kDa polypeptide associated with beta2-microglobulin expressed on cortical thymocytes (strongly), Langerhans cells, dendritic cells and some T cell leukemias and lymphomas. |                  |
| Other names           | T6, Leu-6, HTA1, FCB6                                                                                                                                                                                                                                                |                  |
| Immunogen             | Human thymocytes                                                                                                                                                                                                                                                     |                  |
| Entrez Gene ID        | 909                                                                                                                                                                                                                                                                  |                  |
| Gene name             | CD1A                                                                                                                                                                                                                                                                 |                  |
| NCBI Full Gene Name   | CD1a molecule                                                                                                                                                                                                                                                        |                  |
| UniProt ID            | P06126                                                                                                                                                                                                                                                               |                  |
| Preparation           | Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.                                                                                        |                  |
| Formulation           | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                                                                                                                                                                              |                  |
| Storage and handling  | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.                                                                                                                                                                                             |                  |
| Images and References | <a href="http://www.exbio.cz">www.exbio.cz</a>                                                                                                                                                                                                                       |                  |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.