

Technical Data Sheet

Product

Anti-ARHGEF4 PE

Cat. Number/Size

1P-723-C025

0.025 mg

1P-723-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	ARHGEF4
Clone	ARHGEF-08
Format	PE
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody ARHGEF-08 recognizes human intracellular protein ARHGEF4 / ASEF1, approx. 80 kDa guanine nucleotide exchange factor specific for Rac1 and Cdc42. The epitope is located in the rough at the region of SH3 domain.
Other names	ASEF1, GEF4, STM6
Immunogen	Recombinant fragment of human ARHGEF4 (amino acids 143-271)
Entrez Gene ID	50649
Gene name	ARHGEF4
NCBI Full Gene Name	Rho guanine nucleotide exchange factor 4
UniProt ID	Q9NR80

Concentration	0.1 mg/ml
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.