



	Technical Data Sheet	
Product	Anti-Hu CD11b Purified Low Endotoxin	
Cat. Number/Size	12-681-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD11b	
Clone	ICRF44	
Format	Purified Low Endotoxin	
Reactivity	Human, Non-human primates	
Application	FC (QC tested), IP, IHC-F, ICC, FUNC	
Application details	Functional application: Blocking of granulocyte activation. Flow cytometry: Recommended dilution: 1-12 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody ICRF44 recognizes an extracellular epitope of CD11b (Mac-1alpha), a 165-170 kDa type 1 transmembrane protein mainly expressed on monocytes, granulocytes and NK-cells.	
Other names	Mac-1, Integrin alpha M, ITGAM, CR3A, MO1A, MAC1A	
Workshop	HLDA IV: WS Code M047	
Immunogen	Rheumatoid synovial cells and human monocytes.	
Entrez Gene ID	3684	
Gene name	ITGAM	
NCBI Full Gene Name	integrin subunit alpha M	
UniProt ID	P11215	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.