

Technical Data Sheet

Product	Anti-Hu CD82 Purified Low Endotoxin
Cat. Number/Size	12-559-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	CD82
Clone	C33
Format	Purified Low Endotoxin
Reactivity	Human
Application	FC (QC tested), IP, WB, IHC(P), ICC, FUNC
Application details	Functional application: In human MOLT-4 T-cell line the antibody C33 inhibits syncytium formation induced by coculture with human T-cell leukemia virus type 1 (HTLV-1)-positive human T-cell lines. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody C33 recognizes an extracellular/luminal epitope of CD82, a widely expressed cell surface protein of the tetraspanin family. CD82 is also found in endosome/lysosome compartments.
Other names	R2, 4F9, C33, IA4, ST6, GR15, KAI1, SAR2, TSPAN27
Immunogen	C91/PL (human HTLV-1+ T cell line)
Entrez Gene ID	3732
Gene name	CD82
NCBI Full Gene Name	CD82 molecule
UniProt ID	P27701

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.