

Technical Data Sheet

Product **Anti-Ms LPAM-1 Purified Azide Free**
Cat. Number/Size **10-747-C100 0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen LPAM-1
Clone DATK32
Format Purified Azide Free
Reactivity Mouse
Application FUNC, FC (QC tested), IP, IHC(F)
Application details Functional application: Blocking of cell adhesion.
Isotype Rat IgG2a kappa
Specificity The rat monoclonal antibody DATK32 recognizes an extracellular epitope of integrin alpha 4 (CD49d) and integrin beta7 (Ly69) components of mouse LPAM-1 complex, which is expressed on the majority of peripheral lymphocytes, as well as on subsets of thymocytes and bone marrow cells.
Other names Integrin alpha 4 beta 7, CD49d / Ly69
Immunogen TK1 cells

Concentration 1 mg/ml
Preparation Purified by protein-G affinity chromatography.
Formulation Phosphate buffered saline (PBS), pH 7.4
Storage and handling Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.