

Technical Data Sheet

Product	Anti-Tb10.3 (M. tuberculosis) Purified	
Cat. Number/Size	11-187-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Tb10.3 (M. tuberculosis)
Clone	PAb (187)
Format	Purified
Reactivity	Mycobacterium tuberculosis
Application	WB, ELISA
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody to Tb10.3 reacts with Mycobacterium tuberculosis protein Tb10.3 (EsxR).
Other names	Rv3019c, EsxR
Immunogen	Recombinant Tb10.3 protein (produced in E. coli)
Entrez Gene ID	888926
Gene name	esxR
NCBI Full Gene Name	ESAT-6 like protein EsxR
UniProt ID	P9WNI9

Concentration	1 mg/ml
Preparation	Purified by ligand affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.