



	Technical Data Sheet	
Product	Anti-Hu CD178 Purified	
Cat. Number/Size	11-754-C025	0.025 mg
	11-754-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD178	
Clone	NOK-1	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB, ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody NOK-1 recognizes an extracellular epitope of CD178 / Fas-L, an approximately 40 kDa transmembrane glycoprotein expressed on neutrophils, monocytes, and activated T and NK cells.	
Other names	APTL, FASL, FAS-ligand, CD95L, TNFSF6, TNLG1A	
Workshop	HLDA VII: 70322	
Immunogen	L5178Y mouse T lymphoma cells expressing recombinant human CD178	
Entrez Gene ID	356	
Gene name	FASLG	
NCBI Full Gene Name	Fas ligand	
UniProt ID	P48023	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.