

Technical Data Sheet

Product	Anti-Hu SIGLEC8 Purified	
Cat. Number/Size	11-846-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	SIGLEC8
Clone	7C9
Format	Purified
Reactivity	Human
Application	FC (QC tested), ELISA, IP
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody 7C9 recognizes an extracellular epitope of SIGLEC8, an eosinophil marker, expressed e.g. in lymph nodes and spleen.
Other names	SAF2
Immunogen	Extracellular domain of human SIGLEC8 fused with Fc fragment of human IgG1
Entrez Gene ID	27181
Gene name	SIGLEC8
NCBI Full Gene Name	sialic acid binding Ig like lectin 8
UniProt ID	Q9NYZ4

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.