



	Technical Data Sheet	
Product	Anti-Hu CD154 Purified Low Endotoxin	
Cat. Number/Size	12-781-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD154	
Clone	24-31	
Format	Purified Low Endotoxin	
Reactivity	Human, Non-human primates	
Application	FC (QC tested), ICC, FUNC	
Application details	Functional application: Blocking. Flow cytometry: Recommended dilution: 2-6 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody 24-31 detects an extracellular epitope of CD154 / CD40L (CD40-ligand), a 39 kDa cell surface type II glycoprotein expressed predominantly on activated CD4+ lymphocytes.	
Other names	CD40L, CD40 ligand, TNFSF5, Ly62, gp39, T-BAM	
Immunogen	human CD154 fusion protein	
Entrez Gene ID	21947	
Gene name	CD40LG	
NCBI Full Gene Name	CD40 ligand	
UniProt ID	P27548	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.