



	Technical Data Sheet	
Product	Anti-Hu CD263 FITC	
Cat. Number/Size	1F-518-C025	0.025 mg
	1F-518-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD263	
Clone	TRAIL-R3-02	
Format	FITC	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG1	
Specificity	The antibody TRAIL-R3-02 reacts with TRAIL-R3, a 35 kDa GPI-anchored extracellular membrane protein expressed mainly on neutrophils.	
Other names	TNFRSF10C, DcR1, TRAIL-R3, TRID, LIT	
Workshop	HLDA IX	
Immunogen	TRAIL-R3-hlgGhc fusion protein	
Entrez Gene ID	8794	
Gene name	TNFRSF10C	
NCBI Full Gene Name	TNF receptor superfamily member 10c	
UniProt ID	O14798	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.