



	Technical Data Sheet	
Product	Anti-FOLR2 Purified	
Cat. Number/Size	11-696-C025	0.025 mg
	11-696-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	FOLR2	
Clone	EM-35	
Format	Purified	
Reactivity	Human	
Negative species	Mouse	
Application	FC (QC tested), IP, WB, ICC	
Application details	Western blotting: Non-reducing conditions. Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody EM-35 recognizes an extracellular epitope on FOLR2, a 30-40 kDa cell surface protein serving as a receptor for folic acid.	
Other names	Folate receptor beta, FBP, FR-P3, FR-BETA, BETA-HFR, FBP/PL-1	
Immunogen	BW5147alpha,beta- cells	
Entrez Gene ID	2350	
Gene name	FOLR2	
NCBI Full Gene Name	folate receptor beta	
UniProt ID	P14207	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.